

## HKO Hinged Knee Orthosis – Anterior Closure

**NOTE:** Instructions cover multiple product families, images may vary.



HKO 11" with two straps

HKO 16" with four straps



Figure 1



Figure 2

**NOT MADE WITH  
natural rubber latex**

Please read the following instructions carefully and completely before use. Correct application is required to ensure proper function of the device.

### Intended for single patient use.



**WARNING:** If your pain increases or persists or you experience any unusual reactions while using this product, consult your healthcare provider immediately.

### APPLICATION INSTRUCTIONS

1. Center brace behind knee and wrap around leg (Figure 1).
2. Ensure patella opening is centered over the kneecap. Secure wrap using hook-and-loop contact closures (Figure 1).
3. Wrap and secure the loop-and-lock leg straps (Figure 2).

**NOTE:** If applying the 16" HKO, fasten straps located above and below the knee first.

4. Adjust loop-and-lock leg straps to achieve desired fit and level of compression.

**NOTE:** Please refer to sew-in label on device for care instructions.

Please visit our website at **[www.birdcronin.com](http://www.birdcronin.com)** for additional product information or contact a customer care representative at **800-328-1095**.

**NOTICE:** Every effort has been made to obtain the maximum compatibility of function, strength, durability and comfort with this product, however, there is no guarantee that injury will be prevented with its use.

