Össur Americas
27051 Towne Centre Drive
Foothill Ranch, CA 92610, USA
Tel: +1 (949) 382 3883
Tel: +1 800 233 6263
Fax: +1 800 831 3160
ossurusa@ossur.com

Össur Canada
2150 – 6900 Graybar Road
Richmond, BC
V6W OAS, Canada
Tel: +1 604 241 8152
Fax: +1 866 441 3880

Össur Europe BV
Ekkersrijt 4106-4114
P.O. Box 120
5690 AC Son en Breugel
The Netherlands
Tel: +800 3539 3668
Tel: +31 499 462840
Fax: +31 499 462841
info-europe@ossur.com

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EN – Caution: This product has been designed and tested as a single use product, and is not intended for multiple use. In addition to increasing the risk of infection, used product components may contain small unnoticeable defects that may lead to non-optimal performance, and even mechanical failure. If any problems occur with the use of this product, immediately contact your medical professional.

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

REBOUND® CARTILAGE INDICATIONS FOR USE

- Articular cartilage defect repair needing unicompartmental load reduction or range of movement restriction.
- Cartilage knee trauma benefiting from unloading or movement restriction and needing pain relief.
- May also be beneficial for other unicompartmental knee conditions that require unloading, movement restriction and pain relief such as:
 - Meniscal repair.
 - Avascular necrosis.
 - Condylar bone marrow lesions (bone bruises)
- Circumference limitation.
 - Thigh circumferences smaller than 13.5" (34.3 cm) or larger than 29" (73.5 cm).
 - Calf circumferences smaller than 10" (25 cm) or larger than 24" (61.5 cm).

CAUTION

- For patients with symptomatic popliteal cyst, peripheral vascular disease, neuropathy, and sensitive skin, close physician supervision is recommended. For patients with sensitive skin, it is recommended to use the Rebound® Cartilage as tolerated. If redness appears decrease the usage of the brace. Increase usage slowly as the skin tolerates the brace.
- The Rebound® Cartilage is field adaptable to accommodate patient comfort issues. Modifications to flexion or extension stops, hinge and shell contour should only be made by the health care professional. The healthcare professional is responsible for advising patient on use of the Cartilage Protection Strap™.

INITIAL APPLICATION AND ADJUSTMENT

1. Place the brace on patient's leg. Align the center of hinge at mid patella and midline of A/P.
2. Place button of lower Quick Fit Buckle into appropriate keyhole above Buckle Stability Shelf (**Figure 2**) and close carefully using the palm of the hand to snap the buckle closed. The button is color matched with an identical blue color in the keyhole for intuitive strap selection.
3. Place button of upper Quick Fit Buckle into appropriate keyhole and lose carefully using the palm of the hand. The button is color

matched with an identical yellow color in the keyhole for intuitive strap selection.

4. Adjust thigh and calf strap tension and trim straps as needed.
5. Have patient bend leg at 45° and set both Cartilage Protection Straps™ to “0”. Peel loose the sleeve hooks and adjust the strap tension starting at the calf (**Figure 3**). Both straps are adjusted correctly when the entire strap has contact with the patient’s leg while not applying any unloading force. If the brace is used with ROM restriction, do not bend the knee beyond the prescribed restriction and adjust the tension of both straps within the prescribed ROM restriction. In case of any questions please contact your physician.

USER APPLICATION

1. Open upper and lower Quick Fit Buckles.
2. Loosen both of the Cartilage Protection Straps™ to the “0” position.
3. Place the brace on patient’s leg. Align the center of hinge at mid patella and midline of A/P.
4. Place button of lower Quick Fit Buckle into appropriate keyhole above Buckle Stability Shelf (**Figure 2**) and close carefully using the palm of the hand to snap the buckle closed.
5. Bend knee to an 80° angle with foot flat on the floor. If the brace is used with ROM restriction, do not bend the knee beyond the prescribed restriction.
6. Place button of upper Quick Fit Buckle into appropriate keyhole and close carefully using the palm of the hand.
7. Turn both dosing dials clockwise to tighten the Cartilage Protection Strap™ to prescribed position (**Figure 4**). Practitioner will need to determine optimum unloading based on patient’s pain relief feedback. The strap tension should be set to provide optimal pain relief once the dosing dial is set to the setting of “5”.
8. Fully extend leg in a sitting position and make sure that brace is still aligned correctly on leg (**Figure 5**). If the brace is used with extension restriction, do not bend the knee beyond the prescribed restriction and adjust the tension of both straps within the prescribed ROM restriction. In case of any questions please contact your physician.
9. If the Cartilage Protection Strap™ has been pulled to its maximum and patient requires more pain relief, both straps may require shortening at the lower and upper Quick Fit Buckles.
10. Tamper proof keyhole locks may be inserted into calf and thigh shells if desired to ensure straps cannot accidentally be disengaged from shells (**Figure 6**).

NOTE: Have patient wear brace for a few minutes to make sure straps are at comfortable tension.

BRACE REMOVAL

1. Sit with the knee flexed to 90°. If the brace is used with ROM restriction, do not bend the knee beyond the prescribed restriction and remove the brace within the prescribed ROM restriction. In case of any questions please contact your physician.
2. Turn both dosing dials counterclockwise to release tension on the Cartilage Protection Strap™ and adjust back to a setting of “0” (**Figure 7**).
3. Open both upper and lower Quick Fit Buckle.
4. Remove Quick Fit Buckles from keyholes.

5. Remove brace from leg.

ONE-TIME STRAP LENGTH ADJUSTMENT (RECOMMENDED)

1. Remove alligator clip and trim strap to desired length. Alligator clips should be positioned near the Quick Fit buckles, away from popliteal region (**Figure 8**). **NOTE:** to comfortably adjust the Cartilage Protection Straps™, first peel back the sleeve hooks so that they do not interfere with the strap loop.
2. Repeat for all straps.

RANGE OF MOTION ADJUSTMENT (ROM)

The device comes with a 0° extension stop installed. Additional extension stops of 5°, 10°, 15° and 30° have been provided. Flexion stops of 0°, 30°, 60°, 90° are also included.

NOTE: Optimal corrective forces are applied at full extension. If patient requires an extension stop change, make sure to re-adjust tension in Cartilage Protection Straps™ if more unloading is required.

TO CHANGE RANGE OF MOTION

1. Peel condyle pad from hinge to expose screw.
2. Remove retaining screw.
3. Slide 0° stop out of hinge.
4. Insert desired extension stop/flexion stop.
5. Secure with retaining screw (**DO NOT OVER TIGHTEN**).

BRACE MAINTENANCE

Rebound® Cartilage is designed for minimum maintenance and care. Follow these easy maintenance steps to keep brace in good condition.

CLEANING

- Liner sleeve and liners may be wiped down while on the brace with disinfectant liquid to prevent build-up of harmful bacteria and wiped with an antibacterial detergent once per week.
- For more thorough cleaning: liners, straps and strap pads may be removed for cleaning.
- Recommended that sleeve hooks are not removed and stay attached to the liner sleeve straps when cleaning.
- Hand-wash using mild non-detergent soap.
- **DO NOT MACHINE WASH.**
- **DO NOT PLACE IN DRYER.**
- Rinse thoroughly and air dry.
- To aid in strap application after cleaning, see **Figure 9** for final brace configuration.

HINGE

- For safety reasons, hinge should only be disassembled by a health care professional.
- Check regularly for foreign materials in hinge (i. e. dirt or grass).
- Clean dirt from hinge using tap water and/or compressed air.

NOTE: If brace is used in salt water or chlorinated water, rinse well with fresh tap water and air dry. A regular service routine will ensure hinge, liners, straps and pads are in optimal working condition.

BRACE ACCESSORIES

Össur offers accessories that meet specific needs of the Rebound® Cartilage. These accessories can be purchased through the brace provider or Össur Customer Service Department:

- 
- Extra Set of Liner Sleeves.
 - Extra Set of Liners.
 - Replacement Straps.
 - Strap pads.
 - Doeskin Calf Liner.
 - Suspension Strap Kit (**NOTE:** Always recommended with a Doeskin calf liner).
 - Undersleeve (**NOTE:** Optimal solution for patients with sensitive skin and are using regular liners on the brace).

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