FORA[®] G20

Blood Glucose Monitoring System

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Owner's Manual
Manual del Propietario

Dear Owner of FORA G20 Blood Glucose Monitoring System:

Thank you for purchasing the **FORA G20** Blood Glucose Monitoring System. This manual provides important information to help you to use the system properly. Before using this product, please read and follow the instructions in this Owner's Manual.

Regular monitoring of your blood glucose levels can help you and your doctor gain better control of your diabetes. Due to its compact size and easy operation, you can use the **FORA G20** Blood Glucose Monitoring System to easily monitor your blood glucose levels.

If you have questions, just call ForaCare Customer Care Service Center toll-free at 1-888-307-8188 (8:30 am - 5:00 pm PST, Monday to Friday). You can also visit www.foracare.com.

TABLE OF CONTENTS

IMPORTANT SAFETY INSTRUCTIONS	1
Read Before Use	1
Important Information	3
Intended Use	4
Test Principle	4
Contents of System	5
Appearance and Key Functions of the Meter	6
Meter Display	7
Test Strip	8
Setting the Meter and Deleting the Memory	9
CHECKING THE SYSTEM WITH CONTROL SOLUTIONS	12
How to Perform a Control Solution Test	12
TESTING YOUR BLOOD	15
The Two Measuring Modes	15
Preparing a Blood Glucose Test	16
Preparing the Puncture Site	16
Alternative Site Testing	17
Performing a Blood Glucose Test	19
VIEWING THE METER MEMORY	21
Viewing Results	21

Viewing Day Average Results	22
VIEWING RESULTS ON A PERSONAL COMPUTER	23
CLEANING AND DISINFECTING YOUR METER	24
Caring for Your Meter	24
Disinfecting Procedures	25
DISPLAY MESSAGES AND PROBLEM-SOLVING GUIDE	27
Troubleshooting	29
DETAILED INFORMATION	30
Reference Values	30
Comparing Meter and Laboratory Results	30
SPECIFICATIONS	32

IMPORTANT SAFETY INSTRUCTIONS

Read Before Use

- The meter and lancing device are for single patient use only. Do not share them with anyone including your family members! Do not use on multiple patients!
- All parts of the kit are considered biohazardous and can potentially transmit infectious diseases, even after you have performed cleaning and disinfection.

For more information, please visit

- “FDA Public Health Notification: Use of Fingerstick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens: Initial Communication” (2010)
<http://www.fda.gov/MedicalDevices/Safety/AlertsandNotices/ucm224025.htm>
 - “CDC Clinical Reminder: Use of Fingerstick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens” (2010)
<http://www.cdc.gov/injectionsafety/Fingerstick-DevicesBGM.html>
1. Use this device **ONLY** for the intended use described in this manual.
 2. Do **NOT** use accessories which are not specified by the manufacturer.
 3. Do **NOT** use the device if it is not working properly or if it is damaged.
 4. Do **NOT** use the equipment in places where aerosol sprays are being used or where oxygen is being administered.

5. Do **NOT** under any circumstances use the device on newborns or infants.
6. This device does **NOT** serve as a cure for any symptoms or diseases. The data measured is for reference only. Always consult your doctor to have the results interpreted.
7. Before using this device to test blood glucose, read all instructions thoroughly and practice the test. Carry out all the quality control checks as directed.
8. Keep the device and testing equipment away from young children. Small items such as the battery cover, batteries, test strips, lancets and vial caps are choking hazards.
9. Use of this instrument in a dry environment, especially if synthetic materials are present (synthetic clothing, carpets etc.) may cause damaging static discharges that may cause erroneous results.
10. Do **NOT** use this instrument in close proximity to sources of strong electromagnetic radiation, as these may interfere with the accurate operation.
11. Proper maintenance and periodically control solution test are essential to the longevity of your device. If you are concerned about your accuracy of measurement, please contact the local customer service or place of purchase for help.

KEEP THESE INSTRUCTIONS IN A SAFE PLACE

Important Information

- Severe dehydration and excessive water loss may cause readings which are lower than actual values. If you believe you are suffering from severe dehydration, consult a healthcare professional immediately.
- If your blood glucose results are lower or higher than usual, and you do not have symptoms of illness, repeat the test. If you have symptoms or continue to get results which are higher or lower than usual, follow the treatment advice of your healthcare professional.
- Use only fresh whole blood samples to test your blood glucose. Using other substances will lead to inaccurate results.
- If you are experiencing symptoms that are inconsistent with your blood glucose test results and you have followed all the instructions given in this owner's manual, contact your healthcare professional.
- We do not recommend using this product on severely hypotensive individuals or patients who are in shock. Readings which are lower than actual values may occur for individuals in a hyperglycemic-hyperosmolar state, with or without ketosis. Please consult your healthcare professional before use.

Intended Use

The **FORA G20** Blood Glucose Monitoring System is intended for use in the quantitative measurement of glucose in fresh capillary whole blood from the finger and the following alternative sites: palm, forearm or upper-arm. It is intended for use by people with diabetes mellitus at home as an aid in monitoring the effectiveness of diabetes control program. It is not intended for the diagnosis of or screening for diabetes mellitus, and is not intended for use on neonates. It is intended to be used by a single person and should not be shared.

The alternative site testing in the **FORA G20** Blood Glucose Monitoring System can be used only during steady-state blood glucose conditions.

This meter has some speaking functions but is not intended for use by the visually impaired.

Test Principle


Your system measures the amount of sugar (glucose) in whole blood. The glucose testing is based on the measurement of electrical current generated by the reaction of glucose with the reagent of the strip. The meter measures the current, calculates the blood glucose level, and displays the result. The strength of the current produced by the reaction depends on the amount of glucose in the blood sample.

NOTE:

- Check your system to ensure that it is unopened prior to use and that it contains all parts listed above. If anything is missing or damaged, please return your system to the place of purchase.

Contents of System

The **FORA G20** system kit includes:

① **Blood Glucose Meter**

⑤ **Daily Log Book**

② **Owner's Manual**

⑥ **CR2032 lithium battery**

③ **Protective Wallet**

④ **Quick Start User Guide /**


Warranty Card

Test strips, control solution, sterile lancets and lancing device may be purchased separately. Please make sure you have those items needed for a blood glucose test beforehand.

NOTE:

- If any items are missing from your kit or opened prior to use, or to obtain test strips and control solutions, please contact the place of purchase or call Customer Service for assistance.

Appearance and Key Functions of the Meter


- ① **Test Slot**
- ② **Display Screen**
- ③ **Main Button**
- ④ **Strip-Ejection Button**
- ⑤ **Data Port**
- ⑥ **C Button (Control Button)**
- ⑦ **SET Button**

NOTE:

- The meter will turn off automatically after 180 seconds of idle time, or you can press and hold the main button for 3 seconds to turn off the meter.


Meter Display

- | | |
|------------------------|--------------------|
| ① Low Battery Symbol | ⑨ Time |
| ② Test Trip Symbol | ⑩ Alarm Indicator |
| ③ Blood Drop Symbol | ⑪ Day Average |
| ④ Test Result Display | ⑫ Memory Symbol |
| ⑤ Error Message Symbol | ⑬ Measurement Unit |
| ⑥ Ketone Warning | |
| ⑦ CTL Symbol | |
| ⑧ Date | |

Test Strip

Your system measures the amount of sugar (glucose) in whole blood. Blood is applied to the absorbent hole on the test strip and is automatically drawn into the reaction cell where the reaction takes place.

ATTENTION:


The front side of test strip should face up when inserting test strip.

Test results might be wrong if the contact bar is not fully inserted into the test slot.


NOTE:

- The FORA **G20** Blood Glucose Monitoring System should only be used with FORA **G20** Test Strips. Using other test strips with this monitoring system can produce inaccurate results.

Setting the Meter and Deleting the Memory



Start with the meter off (no test strip inserted). Press **SET**.


1. Setting the date

With the year flashing, press Main button until the correct year appears. Press **SET**.

With the month flashing, press Main button until the correct month appears. Press **SET**.






With the day flashing, press Main button until the correct day appears. Press **SET**.

2. Setting the time format


Press Main button to select the desired time format — 12h or 24h. Press **SET**.

3. Setting the hour and minute

With the hour flashing, press Main button until the correct hour appears. Press **SET**. With the minute flashing, press Main button until the correct minute appears. Press **SET**.



4. Deleting the memory

With “dEL” and a flashing “**M**” on the display, press Main button and select “no” to keep the results in memory then press **SET** to skip.


To delete all the results, press Main button and “yes” and “**M**” are displayed on the meter, press **SET** to delete the memory.

If you do not want to delete the results, do not press the Main button but press the **SET** button again to turn off the meter.

5. Setting the reminder alarm

You may set up any or all of the reminder alarms (1-4). When the meter displays “On” or “OFF” and “**M**”, press Main button to turn on or turn off the first reminder alarm. Press **SET** to confirm.

When the alarm is enabled with hour flashing, press Main button to increase an hour, followed by pressing **SET** to confirm and go to minute setting.

Press **SET** to confirm and go to the next alarm setting. If you do not want to set an alarm, press **SET** to skip this step.

If you want to turn off an alarm, find the alarm number by pressing **SET** in the setting mode, and then press Main button to change from “ON” to “OFF”.

Congratulations! You have completed all settings!

NOTE:

- These parameters can **ONLY be changed** in the setting mode.
- If the meter is idle for 3 minutes during the setting mode, it will switch off automatically.

CHECKING THE SYSTEM WITH CONTROL SOLUTIONS

Our control solutions contain a known amount of glucose that react with test strips. By comparing the result of your control solution test with the expected range printed on the test strip vial label, you can check if the meter and test strips are working properly, and if you're performing the test correctly. It is very important that you perform this simple check routinely to make sure you get accurate results.

Do a control solution test when:

1. you first receive the meter,
2. at least once a week to routinely check the meter and test strips,
3. you begin using a new vial of test strips,
4. you suspect the meter or test strips are not working properly,
5. your blood glucose test results are not consistent with how you feel, or if you think the results are not accurate,
6. practicing the testing process, or
7. you dropped or think you have damaged the meter.

How to Perform a Control Solution Test





Test strips, control solutions, or sterile lancets may not be included in the kit (please check the contents on your product box). They can be purchased separately. Please make sure you have all items prior to testing your blood glucose. Please contact the Customer Care Line at 1-888-307-8188 or the place of purchase for availability.

1. Insert test strip

Insert a test strip into the test slot with the contact bars end first and facing up. (Contact bars must be inserted all the way into the meter or you may get an inaccurate test result.) The meter turns on automatically and displays the following in sequence:

When the “” appears on the display, press the Main button and “CTL” will appear on the display. When the “CTL” sign is displayed, the meter will not store your test result in the memory. If you decide not to perform a control solution test, press the Main button again and the “CTL” sign will disappear.

CAUTION:

- Every time you perform a control solution test you must enter “QC” test mode so that the test result will not be stored in the meter memory.

2. Apply control solution

Shake the control solution vial thoroughly before use. Squeeze out a drop and wipe it off, and then squeeze out another drop and place it on the tip of the vial cap. Hold the meter in an angled position so that the absorbent hole of the test strip touches the drop. Once the confirmation window fills completely, the meter will begin counting down. To avoid contaminating the control solution, do not directly apply control solution onto the strip.

3. Read and compare the results

After the meter counts to 0, the result of the control solution test will appear. Compare the result with the range printed on the test strip vial. The result should fall within this range.

Out-of-range results

If you continue to have test results fall outside the range printed on the test strip vial, the meter and strips may not be working properly. Do **NOT** test your blood. Contact customer service or place of purchase for help.

TESTING YOUR BLOOD





Before testing, be sure to read this section and the test strip insert found in the test strip box carefully. Make sure you have all the items needed for the test.

The Two Measuring Modes

The meter provides you with two modes for measuring: General and QC.

MODES	USE WHEN
General tests (not display)	any time of day without regard to time since last meal
QC (CTL)	testing with the control solution

You can switch between the two modes by:

1. Start with the meter switched off. Insert a test strip to turn on the meter. The screen will display a “”.
2. Press Main button to switch between General and QC modes.

Preparing a Blood Glucose Test

Sharing the lancing device and lancets may increase the risk of contracting infectious diseases. The lancing device must not be used for more than one person.

To reduce the chance of infection:

1. The meter and lancing device are for single use. Never share a lancet or the lancing device.
2. Avoid getting hand lotion, oils, dirt, or debris in or on the lancets and the lancing device.
3. Wash and dry your hands thoroughly after handling the meter, lancing device and test strips to prevent infection. For more information, please refer to the “Cleaning and Disinfection” section.
4. If the meter is being operated by a second person who is providing testing assistance to the user, the meter and lancing device should be decontaminated prior to use by the second person.

Preparing the Puncture Site

Rubbing the puncture site before blood extraction has a significant influence on the accuracy of your test results.

When the puncture site was rubbed prior to blood extraction, the test result will be more precise.


Please follow the suggestions below before obtaining a drop of blood:

- **Wash and dry your hands before starting.**
- Select the puncture site on fingertips.
- Rub the puncture site for about 20 seconds before penetration.
- Clean the puncture site using cotton moistened with 70% alcohol and **let it air dry.**

- **Fingertip testing**

Hold the lancing device firmly against the lower side of your fingertip. Press the release button to puncture your finger. You will hear a click indicating that the puncture is complete.


- **Blood from sites other than the fingertip**

Replace the lancing device cap with the clear cap for AST. Pull the cocking control back until it clicks. When lancing the forearm, upper arm or hand, avoid lancing the areas with obvious veins because of excessive bleeding.

NOTE:

- Choose a different spot each time you test. Repeated punctures at the same spot may cause soreness and calluses.
- It is recommended that you discard the first drop of blood as it might contain tissue fluid, which may affect the test result.

Alternative Site Testing

Important:

There are limitations with AST (Alternative Site Testing).

- **AST results should not be used for CGM calibration.**
- **AST results should not be used for insulin dosing calculations.**

Please consult your health care professional before you perform AST.

What is AST?

Alternative site testing (AST) means that people use parts of the body other than the fingertips to check their blood glucose levels. This system provides you to test on the palm, the forearm or the upper arm with the equivalent results to fingertip testing.

What is the advantage?

Your fingertips feel more pain because they are full of nerve endings (receptors). At other body sites, where nerve endings are not so condensed, you will not feel as much pain as you would pricking your fingertip.

When to use AST?

Food, medication, illness, stress and exercise can affect blood glucose levels. Capillary blood at the fingertip reflects these changes faster than capillary blood at other sites. Thus, when testing blood glucose during or immediately after a meal, physical exercise, or any other event, take a blood sample from your finger only. ²⁴

We strongly recommend that you perform AST **ONLY** at the following times:

- In a pre-meal or fasting state (more than 2 hours since the last meal).
- Two hours or more after taking insulin.
- Two hours or more after exercise.

Do **NOT** use AST if:

- You think your blood glucose is low.
- You are unaware of hypoglycemia.
- You are testing for hyperglycemia.
- Your AST results do not match the way you feel.
- Your routine glucose results often fluctuate.


Performing a Blood Glucose Test

1. Insert the test strip to turn on the meter

Wait for the meter to display “”.

2. Select the appropriate measuring mode by pressing Main button

For selecting the measurement mode, please refer to "THE TWO MEASURING MODES".

3. Apply blood sample

Obtain a drop of blood of at least 0.7 μL with a lancing device. When the “” symbol appears on the display, touch the drop of blood with the absorbent hole of the test strip. Do not rub the test strip with your finger or try to apply a smeared sample.

Gently apply the drop of blood to the absorbent hole of the test strip at a tilted angle. Confirmation window should be completely filled if enough blood sample has been applied. Do **NOT** remove your finger until you hear a beep sound.

If the confirmation window is not filled completely before the meter begins to count down, do not add more blood to the test strip. Discard the test strip and start again. If you have trouble filling the test strip, please call your local customer service number for assistance.

NOTE:

- If you do not apply a blood sample within 3 minutes, the meter will automatically turn off. You must remove the test strip and insert it back into the meter to restart the test.

4. Get result

The result of your blood glucose test will appear after the meter counts down to 0. The blood glucose result will be stored in the memory automatically.


WARNING:

- Please do not change your treatment based on the result without first consulting your health care professional.
- Turn the meter off by removing the test strip. Discard the used test strip and lancet properly.
- Wash your hands thoroughly with soap and water after handling the meter, lancet device and test strips to avoid contamination.

VIEWING THE METER MEMORY

Your Meter stores the 450 most recent blood glucose test results with date and time in its memory. You can review the test results with these easy steps.


Viewing Results

1. Enter the memory mode


With the meter turned off, press the Main button. The first test result will appear, indicating that you are in the memory mode.

When using the meter for the first time or if the results have been deleted, “- - -” will appear, indicating that there are no test results in the memory.

2. Recall test results

After the last test result, the most recent test result appears with the date and time. Press the Main button to review your last 450 test results in order. When the memory is full, the oldest result will be deleted when the newest is added.


3. Exit the memory mode

Press the Main button for three seconds to turn off the meter or after the last result, the display will show “End” and it will turn off automatically.

NOTE:


- The results of control solution tests are **NOT** stored in the memory. Stored results are blood glucose results only.

Viewing Day Average Results





1. Press and release Main button

Keep pressing Main button for 3 seconds until the flashing "DAY AVG" appears. Release Main button and then your 7-day average result measured in general mode will appear on the display.

2. Press Main button to review

14-, 21-, 28-, 60- and 90- day average results stored in each measuring mode.

3. Exit the meter memory

Keep pressing the Main button and the meter will turn off after displaying the last test result.

VIEWING RESULTS ON A PERSONAL COMPUTER

Results in memory can be transmitted to a personal computer. Health Care Software System and an interface cable are required before installation. The software can be downloaded from the home page of ForaCare, Inc. An interface cable is required for data transmission, but is not included in the standard kit. To learn more about Health Care Software System or to obtain an interface cable separately, please contact Customer Care Line at 1-888-307-8188 for help.

1. Install Software

Install Health Care Software System on your Computer by following the instructions provided on ForaCare's website, www.foracare.com

2. Connect to Personal Computer

Connect the interface cable to a serial port of your computer. With the meter turned off, connect the internal cable to the data port of the meter. Then "PC" will appear on the display, indicating that the meter is ready to transmit data.

3. Transmit Data

Follow the instructions to transmit data from the meters memory to your computer. Results transmitted will include date and time. Remove the cable and the meter will automatically turn off.

NOTE:

- While the meter is connected to the PC, it is unable to perform a blood glucose test.

CLEANING AND DISINFECTING YOUR METER

Caring for Your Meter

To avoid the meter and test strips attracting dirt, dust or other contaminants, please wash and dry your hands thoroughly before use.

Why the cleaning and disinfection should be performed

Cleaning and disinfection are different. Cleaning is the process of removing dirt (e.g. food debris, grease, dust), disinfection is the process of killing germs (e.g. bacteria and viruses).

When to clean and disinfect the meter

Clean the meter when you see any dirt on it. You should disinfect the meter at least once a week to prevent infection.

How to clean and disinfect the meter

The meter must be cleaned prior to the disinfection. Use one disinfecting wipe to clean exposed surfaces of the meter thoroughly and remove any visible dirt or blood or any other body fluid with the wipe. Use a second wipe to disinfect the meter. Do **NOT** use organic solvents to clean the meter.


We recommend for meter cleaning and disinfection you should use the disinfecting wipe/towelette from below. The following product has been shown to be safe for use with the FORA G20 glucose meter.

▶ Micro-Kill+™ (Micro-Kill Plus™) by Medline
(EPA Reg. No. 59894-10-37549)

To obtain disinfecting wipes and other information, please contact Medline at 1-800-MEDLINE (1-800-633-5463) or visit www.medline.com.

Disinfecting Procedures

- 1.** Take out one disinfecting wipe from the package and squeeze out any excess liquid in order to prevent damage to the meter.
- 2.** Wipe all meter's exterior surface display and buttons. Hold the meter with the test strip slot pointing down and wipe the area around the test slot but be careful not to allow excess liquid to get inside. Keep the meter surface wet with disinfection solution for a minimum of 2 minutes for Micro-Kill+™ wipes. Follow the instructions on the package label of disinfecting wipe.

- 3.** Allow the meter surface to dry completely.
- 4.** Discard the used wipes and never reuse them. Wash your hands thoroughly with soap and water after handling the meter, lancing device and test strips to avoid contamination.

Improper system cleaning and disinfection may result in meter malfunction. If you have a question, please contact customer service at 1-888-307-8188 for assistance.

This device has been validated to withstand 260 cleaning and disinfection cycles using the recommended disinfecting wipe/towelette. The tested number of cycles is estimated by 1 cleaning and disinfection cycle per week for 5 years, the expected life of device.

Stop using the meter if you see any signs of deterioration. For example:

- meter can not be turned on,
- LCD display cracks or becomes cloudy,
- buttons no longer function,
- meter outer casing cracks,
- data can not be transmitted to App,
- color or paint/printing on housing is abnormal, or
- scratches or abrasions on meter are higher than acceptable.


Please contact the customer service for a replacement meter if any of the signs of deterioration are noticed.

NOTE:

- Do **NOT** clean and disinfect the meter while performing tests.
- If the meter is being operated by a second person, the meter and lancing device should be decontaminated prior to use by the second person.
- Do **NOT** allow cleaning and disinfecting solution to get in the test slot.
- If you do get moisture in the test strip slot, wipe it away with a corner of tissue.
- Always dry the meter thoroughly before using it.
- Do not spray the meter directly with cleaning solutions especially those containing water (i.e. soapy water), as this could cause the solution to enter the case inside and damage the electronic components or circuitry.

DISPLAY MESSAGES AND PROBLEM-SOLVING GUIDE

The following is a summary of display messages. If your meter displays an error message, please follow the actions for the error message as described in the table below. If the problem persists, please contact your local customer service agent for help.

MESSAGE	WHAT IT MEANS	ACTION
Lo	Appears when glucose < 20 mg/dL (1.1 mmol/L).	This indicates hypoglycemia (low blood glucose.) Please seek immediate medical assistance.
Hi	Appears when glucose > 600 mg/dL (33.3 mmol/L).	This indicates severe hyperglycemia (high blood glucose). Please seek immediate medical assistance.
	The “  ” symbol appears with this message. This means that the battery does not have enough power remaining for a test.	Recharge the battery immediately.
	A used test strip is inserted, or there is an electrical problem.	Repeat test with a new test strip. If the error message re-appears, please contact your local customer service agent for help.
	You may have removed the strip after applying blood to the absorbent hole.	Keep the test strip in the meter after it has absorbed the blood until the meter displays the test result.

	<p>The temperature of the environment, meter, or test strip is outside the system's operating range. You cannot perform a test until the meter and test strip are within the operating range of 10°C to 40°C (50°F to 104°F).</p>	<p>Repeat the test after the meter and test strip are within the operating temperatures.</p>
   	<p>Problem with the meter.</p>	<p>Review the instructions and try again with a new test strip. If the problem persists, please contact your local customer service agent for help.</p>

For questionable results, please see test strip package insert for important information.

Troubleshooting

1. If the meter does not display a message after inserting a test strip:

POSSIBLE CAUSE	WHAT TO DO
Battery exhausted.	Recharge the battery.
Test strip inserted upside down or incompletely.	Insert the test strip with contact bars end first and facing up.
Defective meter or test strips.	Please contact customer service.

2. If the test does not start counting down after applying the sample:

POSSIBLE CAUSE	WHAT TO DO
Defective test strip.	Repeat the test with a new test strip.
Sample applied after automatic switch-off (3 minutes after last user action).	Repeat the test with a new test strip. Apply sample only when flashing “  ” appears on the display.
Defective meter.	Please contact customer service.

3. If the control solution testing result is out of range:

POSSIBLE CAUSE	WHAT TO DO
Error in performing the test.	Read instructions thoroughly and repeat the test again.
Control solution vial was poorly shaken.	Shake the control solution vigorously and repeat the test again.
Expired or contaminated control solution.	Check the expiry date of the control solution.
Control solution that is too warm or too cold.	Control solution, meter, and test strips should be at room temperature (20°C to 25°C / 68°F to 77°F) before testing.
Defective test strip.	Repeat the test with a new test strip.
Meter malfunction.	Please contact customer service.

DETAILED INFORMATION

Reference Values

Blood Glucose

Blood glucose monitoring plays an important role in diabetes control. A long-term study showed that maintaining **blood glucose levels close to normal** can reduce the risk of diabetes complications by up to 60%*¹. The results provided by this system can help you and your healthcare professional monitor and adjust your treatment plan to gain better control of your diabetes.

Time of day	Normal plasma glucose range for people without diabetes (mg/dL)
Fasting and before meal	< 100 mg/dL (5.6 mmol/L)
2 hours after meals	< 140 mg/dL (7.8 mmol/L)

Source: American Diabetes Association (2014). Clinical Practice Recommendations. Diabetes Care, 37 (Supplement 1): S16.

Please work with your doctor to determine a target range that works best for you.

References:

*1: American Diabetes Association position statement on the Diabetes Control and Complications Trial (1993).

Comparing Meter and Laboratory Results

The meter provides you with plasma equivalent results. The result you obtain from your meter may differ somewhat from your laboratory result due to normal variation. Meter results can be affected by factors and conditions that do not affect laboratory results in the same way. To make an accurate comparison between meter and laboratory results, follow the guidelines below.

Before going to the lab:

- Perform a control solution test to make sure that the meter is working properly.
- Fast for at least eight hours before doing comparison tests, if possible.
- Take your meter with you to the lab.

While staying at the lab:

Make sure that the samples for both tests are taken and tested within 15 minutes of each other.

- Wash your hands before obtaining a blood sample.
- Never use your meter with blood that has been collected in a gray-top test tube.
- Use fresh capillary blood only.

You may still have a variation from the result because blood glucose levels can change significantly over short periods of time, especially if you have recently eaten, exercised, taken medication, or experienced stress*². In addition, if you have eaten recently, the blood glucose level from a finger prick can be up to 70 mg/dL (3.9 mmol/L) higher than blood drawn from a vein (venous sample) used for a lab test*³.

Therefore, it is best to fast for eight hours before doing comparison tests. Factors such as the amount of red blood cells in the blood (a high or low hematocrit) or the loss of body fluid (dehydration) may also cause a meter result to be different from a laboratory result.

References:

*2: Surwit, R.S., and Feinglos, M.N.: *Diabetes Forecast* (1988), April, 49-51.

*3: Sacks, D.B.: "Carbohydrates." Burtis, C.A., and Ashwood, E.R. (ed.), *Tietz Textbook of Clinical Chemistry*. Philadelphia: W.B. Saunders Company (1994), 959.

SPECIFICATIONS

Dimension: 87 (L) x 51 (W) x 15 (H) mm

Weight: 42 g (without battery)

Power Source: CR2032 lithium battery

Display: LCD

Memory: 450 measurement results with date and time

External Output: RS-232

Automatic detection of electrode insertion

Automatic reaction time count-down

Temperature warning

Auto turn-off after 3 minutes without action

Operating Condition:

50 °F to 104 °F (10 °C to 40 °C), below 85% R.H.

Storage/Transportation Conditions:

-4 °F to 140 °F (-20 °C to 60 °C), below 95% R.H.

Measurement Units: fixed mg/dL

Measurement Range: 20 to 600 mg/dL (1.1 to 33.3 mmol/L)

This device has been tested to meet the electrical and safety requirements of: IEC/EN 61010-1, IEC/EN 61010-2-101, EN 61326-1, EN 61326-2-6.

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Distributed by Foracare, Inc. / Distribuido por Foracare, Inc.