
# FORA<sup>®</sup> V12 No-Coding

Blood Glucose Monitoring System

**Owner's Manual**

Thank you for choosing the

# **FORA V12 NO-CODING BLOOD GLUCOSE MONITORING SYSTEM**

Dear FORA V12 No-Coding System Owner:

This manual contains important information you must know about the system. Please read it thoroughly and carefully.

For other questions regarding this system, please contact your local Customer Care Line. For all other questions, contact your healthcare professional for assistance.

# IMPORTANT SAFETY INSTRUCTIONS

## READ THIS BEFORE USING

- The meter and lancing device are for single patient use. Do not share them with anyone including other family members! Do not use on multiple patients!
- All parts of the kit are considered biohazardous and can potentially transmit infectious diseases, even after you have performed cleaning and disinfection.

For more information, please visit

1. "FDA Public Health Notification: Use of Fingerstick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens: Initial Communication" (2010)  
<http://www.fda.gov/MedicalDevices/Safety/AlertsandNotices/ucm224025.htm>
2. "CDC Clinical Reminder: Use of Fingerstick Devices on More than One Person Poses Risk for Transmitting Bloodborne Pathogens" (2010) <http://www.cdc.gov/injectionsafety/Fingerstick-DevicesBGM.html>

The following basic safety precautions should always be taken.

1. Close supervision is necessary when the device is used by, on, or near children, handicapped persons or invalids.
2. Use the device only for the intended use described in this manual.
3. Do NOT use accessories which are not supplied by the manufacturer.
4. Do NOT under any circumstances use the device on neonates or infants.

5. Do NOT use the equipment where aerosol sprays are being used, or where oxygen is being administered.
6. Do NOT use device if it is not working properly, or if it has suffered any damage.
7. Read all instructions thoroughly and practice the test before using the product to test your blood glucose. Do all quality control checks as directed and consult with a diabetes healthcare professional.
8. This device does NOT serve as a cure for any symptoms or diseases. The data measured is for reference only.
9. Keep the device and testing equipment away from young children. Small items such as the battery cover, batteries, test strips, lancets and vial caps are choking hazards.
10. Use of this device in a dry environment, especially if synthetic materials are present (synthetic clothing, carpets etc.) may cause damaging static discharges that may cause erroneous results.
11. Do not use this device in close proximity to sources of strong electromagnetic radiation, as these may interfere with the accurate operation.

**KEEP THESE INSTRUCTIONS IN A SAFE PLACE**

# TABLE OF CONTENTS

<b>IMPORTANT INFORMATION</b>	<b>08</b>
<b>ABOUT ALTERNATIVE SITE TESTING (AST)</b>	<b>09</b>
<b>INTRODUCTION TO THE SYSTEM</b>	<b>11</b>
Intended Use	11
Principle of Measurement	11
Contents of the System	12
Appearance and Key Function of the Meter	13
Meter Display	14
Speaking Function	15
Test Strip	16
<b>BEFORE USE</b>	<b>17</b>
Setting the Meter and Deleting the Memory	17
<b>THE TWO MEASURING MODES</b>	<b>20</b>
<b>CONTROL TEST</b>	<b>21</b>
Important Control Solution Test Information	22
Performing a Control Solution Test	23
<b>TESTING BLOOD SUGAR</b>	<b>25</b>
Test Procedure	26
Expected Test Results	31
<b>COMPARING METER AND LABORATORY RESULTS</b>	<b>32</b>
<b>USING THE METER MEMORY</b>	<b>33</b>
Viewing Results On the Meter	33
Viewing Results On a Personal Computer	36

<b>MAINTENANCE</b>	<b>37</b>
Battery Replacement	37
Caring for Your Meter	39
Caring for Your Lancing Device	41
Caring for Your Test Strips	42
<b>PROBLEM-SOLVING GUIDE</b>	<b>43</b>
Special Message	44
Error Message	45
Operational Problem	46
<b>SPECIFICATIONS</b>	<b>48</b>
<b>SUMMARY OF OPERATION</b>	<b>49</b>

## IMPORTANT INFORMATION

- Severe dehydration and excessive water loss may cause readings which are lower than actual values. If you believe you are suffering from severe dehydration, consult a healthcare professional immediately.
- If your blood glucose results are lower or higher than usual, and you do not have symptoms of illness, first repeat the test. If you have symptoms or continue to get results higher or lower than usual, follow the treatment advice of your healthcare professional.
- Use only capillary whole blood sample to test your blood glucose. Using other substances will cause wrong results.
- If you are experiencing symptoms that are inconsistent with your blood glucose test results and you have followed all instructions described in this owner's manual, call your healthcare professional.
- Inaccurate results may occur in severely hypotensive individuals or patients in shock. Readings which are lower than actual values may occur for individuals experiencing a hyperglycemic-hyperosmolar state, with or without ketosis. Critically ill patients should not be tested with blood glucose meters.
- Please refer to your test strip package insert for additional important information.

## Alternative Site Testing (AST)

**Important:** There are limitations for doing AST (Alternative Site Testing).


\* AST results should not be used for CGM calibration.

\* AST results should not be used for insulin dosing calculations.

Please consult your health care professional before you do AST.

### What is AST?

Alternative site testing (AST) refers to when individuals check their blood glucose levels using areas of the body other than the fingertip. This system allows AST at the palm, the forearm, the upper arm, the calf or the thigh with results as reliable as those obtained from fingertip testing.

### What's the advantage?

Fingertips feel pain more readily because they are full of nerve endings (receptors). Since nerve endings are not so condensed at other body sites, tests performed there will be less painful.

### When to use AST?

Food, medication, illness, stress and exercise can affect blood glucose levels.

Capillary blood at the fingertip reflects these changes faster than capillary blood at other sites. Therefore, when testing blood glucose during or immediately after a meal or physical exercise, results from AST and fingertip may be significantly different.

We strongly recommend you do AST in the following intervals:

- In a pre-meal or fasting state (more than 2 hours since the last meal).
- Two hours or more after taking insulin.
- Two hours or more after exercise.

Do **NOT** use AST if:

- You think your blood glucose is low.
- You are unaware of hypoglycemia.
- You are testing for hyperglycemia.
- Your AST results do not correspond with the way you feel.
- Your routine glucose results are often fluctuating.
- You are pregnant.

### **How to increase the accuracy?**

Stimulating blood perfusion by rubbing the puncture site prior to blood extraction has a significant influence on the glucose value obtained. Blood from a site that has not been rubbed exhibits a measurably different glucose concentration than blood from the finger. When the puncture site is rubbed prior to blood extraction, the difference is markedly reduced.

**Please follow suggestions below before puncturing the skin:**

- Rub the puncture site approximately 20 seconds before penetration.
- Use a clear cap (included in the kit) while setting the lancing device.

# INTRODUCTION TO THE SYSTEM

## Intended Use

---

The system is intended for use outside the body (in vitro diagnostic use). It should be used only for testing glucose (sugar) and only with fresh capillary whole blood samples (from the finger, the palm, the forearm, the upper arm, the calf and the thigh). It should not be used for the diagnosis of diabetes or for the testing of newborns. It is intended to be used by a single person and should not be shared.

AST in this system can be used only during steady-state blood glucose conditions described in the section of “About Alternative Site Testing (AST)”.

## Principle of Measurement

---

The test is based on the measurement of electrical current generated by the reaction of glucose with the reagent of the strip. The meter measures the current and displays the corresponding blood glucose level. The strength of the current produced by the reaction depends on the amount of glucose in the blood sample.


This system uses enzyme of GOD: glucose oxidase (GOD) from *A. niger*.

## Contents of the System

---

The system consists of three main products: the blood glucose meter, test strips, and a control solution. These products have been designed, tested, and proven to work together as a system to produce accurate blood glucose test results. Use only the FORA V12 No-Coding Blood Glucose Test Strips and a control solution with the blood glucose meter.

Your system includes:


### PLEASE NOTE

1. Please make sure that all products listed above are contained and sealed in the package before using this system. If you find any imperfection in our products, please return the whole system to the place of purchase.
2. Test strips, control solutions and lancets are optional. They are not included in the standard kit. Please contact the Customer Care Line at 1-888-307-8188 or place of purchase for availability.


## Appearance and Key Function of the Meter

---

- 1 **TEST SLOT** is where you insert the test strip. The meter will turn on automatically after insertion.
- 2 **LCD SCREEN** guides you through the test using symbols and simple messages.
- 3 **MAIN BUTTON** located in front of the meter with "M" on it, is used to turn on the meter, enter the memory mode and control steps of setting.
- 4 **SET BUTTON** (in the back) located in the battery compartment, is used to set up the meter.
- 5 **DATA PORT** located at the side, is for cable connection.

## Meter Display

- 1 TEST STRIP SYMBOL**  
Appears when the meter is turned on.
- 2 BLOOD DROP SYMBOL**  
Flashes when it is ready to apply the sample.
- 3 TEMPERATURE SYMBOL**  
Appears when ambient temperature is outside operating temperature.
- 4 CTL SYMBOL**  
Appears when doing a control test and indicates that the result won't be stored in the memory.
- 5 DATE**
- 6 TIME**
- 7 VOICE SYMBOL**  
Indicates speaking function is on/off.
- 8 LOW BATTERY SYMBOL**  
Appears when the battery power is low.
- 9 TEST RESULT AREA**  
Displays glucose results.
- 10 UNIT OF MEASUREMENT**
- 11 KETONE WARNING**  
Appears when the test result is equal or higher than 240mg/dL.
- 12 DAY AVERAGE**  
Indicates that the displayed test result is an average.
- 13 FACE SYMBOL**  
Appears when memory deletion completed.
- 14 MEMORY SYMBOL**  
Appears when you review the memory.

## Speaking Function

---

When a voice symbol  is on the screen, it means that the meter provides speaking function. It “speaks” step by step that could lead you to test your blood glucose in order. The following table tells you when and what the meter “speaks”. Those messages can also be found in the final section of this manual, “SUMMARY OF OPERATION”.

<b>WHEN</b> does the meter speak?	<b>WHAT</b> does the meter speak?
When you turn on the meter.	Thank you for using this product. Please relax during measurement.
When room temperature is outside operating range, which is 50°F to 104°F (10°C to 40°C).	Room temperature out of range, unable to measure.
When the meter is ready to test. (symbol appears on display).	Please apply blood into the strip.
When the test is completed. (result appears on display).	The blood glucose is (number) milligram per deciliter.
When the test result is outside measurement range, which is 20-600 mg/dL.	The blood glucose out of range, unable to measure.

### NOTE

FORA V12 No-Coding blood glucose monitoring system contains some speaking functions but has not been validated for use by visually impaired users.

## Test Strip

---

This system measures the amount of sugar (glucose) in whole blood. Blood is applied to the absorbent hole of the test strip and is automatically drawn into the reaction cell where the reaction takes place.

The test strip consists of the following parts:

**Contact Bars** .....

Insert this end of the test strip into the meter.  
Push it in firmly until it will go no further.

**Test Strip Handle** .....


Hold this part to insert the test strip into the slot.

**Confirmation Window** .....

This is where you confirm if enough blood has been applied to the absorbent hole of the strip.

**Absorbent Hole** .....

Apply a drop of blood here,  
The blood will be sucked automatically.

Please see section, "Testing Your Blood Sugar", for complete instructions.

## BEFORE USE

### Setting the Meter and Deleting the Memory

Your meter comes with the time, date, memory deletion, voice volume and language preset. If you need to set these parameters, please follow below steps.

**Start with the meter off. Then press the set button located in the battery compartment. The meter is now in the setting mode.**

#### **STEP 1. Set the Year**

After pressing the set button, the year, a blinking number, will be shown on the screen firstly. Press and release the M button to plus one year. You can also keep pushing down the M button to proceed faster. When the correct year is displayed on the screen, press the set button and then a flashing number, which stands for the month, is shown.

#### **STEP 2. Set the Month**

Press and release the M button until you see the correct month. To move faster, keep pushing the M button down. When the desired month is displayed, press the set button and the date will flash.

### STEP 3. Set the Date

Press and release the M button until you see the correct date. To move faster, keep pushing the M button down. When the desired date is displayed, press the set button and then the hour will be shown on the screen.

### STEP 4. Set the Hour

Press and release the M button until you see the correct hour on the screen. To move faster, keep pushing the M button down. When the desired hour is displayed on the screen, press the set button and then the minute will flash.


### STEP 5. Set the Minutes

Press and release the M button to advance one minute. To move faster, hold the M button down. With the correct minute on the display, press the set button and the memory deletion will start flashing.

### STEP 6. Delete Memory



"dEL" and flashing "M" symbol appear on the display. If you don't want to delete memory, press the set button again to skip this step. If you'd like to delete ALL memory, press M button. Both "dEL" and "M" will flash. Press the M button again to delete all memory. The meter then displays "☺", which means that the memory is deleted.

### STEP 7. Select Speaking Function

The meter displays “VOL”, “” and flashing number. Press M button to select speaking volume from 0 to 7. Then press the set button to proceed to the next step.

Number 0 indicates that the speaking function is turned off. “” will not display on LCD during testing. Numbers 1 to 7 indicates speaking volume from low to high. At this time “” will be displayed on LCD during test.

---

### STEP 8. Select Language

L1/L2 and “” display on LCD. Press the M button to select one of the two languages.

Set-up is completed. Press the set button to turn off the meter. “OFF” will be displayed before shut down.

### NOTE

- The time, date, memory deletion, voice volume and languages can ONLY be changed in the setting mode. Therefore, when you perform a glucose testing, those parameters are not possible to be changed.
- Your meter displays 7-, 14-, 21-, 28-, 60- and 90-day averages which you can access from the meter memory. These averages are calculated from results obtained during the 7-, 14-, 21-, 28-, 60- and 90-day preceding the current date and time settings. When the date and time are changed, the 7-, 14-, 21-, 28-, 60- and 90-day averages may change.
- While the meter is in the setting mode, it will turn off automatically without any action in three minutes.


## THE TWO MEASURING MODES

The meter provides you with two modes for measuring, General and CTL.

MODES	USE WHEN
General (Not display)	any time of day without regard to time since last meal
CTL	testing with the control solution

You can switch between each mode by:

1. Start with the meter switched off. Insert a test strip to turn on the meter, the screen will display a flashing blood drop with date and time.
2. Press M to switch to CTL mode.

## CONTROL TEST

FORA control solutions contain a known amount of glucose that reacts with test strips. By comparing your control solution test results with the expected range printed on the test strip vial label, it is able to check that the meter and the test strips are working together as a system and that you are performing the test correctly. It is very important that you do this simple check routinely to make sure you get accurate results.

### **How often should the control solution test be performed?**

- When you use this system to test your blood for the first time, practice the procedure using control solution. When you can do three tests in a row that are within the expected range, you are ready to test your blood.
- For routinely check the meter and test strips, perform a single test for each level of control solution at least once a week.

### **When should the control solution test be performed?**

- When you first get your meter.
- When you begin using a new vial of test strips.
- Whenever you suspect that the meter or test strips are not working properly.
- When your blood glucose test results are not consistent with how you feel, or when you think your results are not accurate.
- When your test strips are exposed to extreme environmental conditions (See Storage section of this manual).
- When you want to practice running the test.
- If you drop the meter.

## Important Control Solution Test Information

---

- Use only FORA control solutions.
- Check the expiration date on the control solution vial. Do not use if expired.
- Control solution, meter, and test strips should come to room temperature (68 to 77°F/20 to 25°C) before testing.
- Shake the vial, discard the first drop of control solution, and wipe off the dispenser tip to ensure a clean sample and an accurate result.
- Use only 90 days after first opening. Record the discard date (date opened plus 90 days) on the control solution vial. Discard after 90 days.
- Store the control solution tightly closed at temperatures 36°F to 86°F (2°C to 30°C). Do not freeze.

### PLEASE NOTE

The control solution range printed on the test strip vial is for control solution only. It is used to test meter and test strip performance. It is not a recommended range for your blood glucose level.

## Performing a Control Solution Test

### STEP 1. Insert Test Strip


Fully insert test strip, face up with control bars end first, into test slot. The meter will turn on automatically and display the followings in sequence:

→ CH" and "  " → a flashing symbol "  " with date and time.

### STEP 2. Press the M Button

While the "  " symbol appears on the display, press the M button so that "CTL" will appear on the display. With the "CTL" sign on the display, the meter will not store your test result in memory. If you decide not to perform a control solution test, press the M button again, and the "CTL" sign will disappear.

### WARNING

- Contact bars must be inserted all the way into the meter or you may receive inaccurate test result.
- Every time you perform a control solution test, you must enter into the "CTL" test mode so that the test result will not be stored in the meter memory. Failure to do so will cause confusion between the blood glucose test result and the control solution test result in system memory.



### STEP 3. Obtain Control Solution

Shake the control solution vial well. Remove the cap from the control solution bottle. Place cap on flat surface. Squeeze the vial, discard the first drop, and wipe off the dispenser tip to prevent contamination. Squeeze the vial again to produce another drop and **place this drop on the top of cap**.

### STEP 4. Apply Control Solution


While holding the meter, move the absorbent hole of the test strip to touch the drop of control solution. Then the drop will be automatically drawn into the test strip. Make sure the confirmation window fills completely. The meter begins to count down.

To avoid contaminating the control solution with the contents of the test strip, you have to place a drop of control solution on a clean surface. Do not directly apply control solution to a strip.

### STEP 5. Read and Compare the Result

After counting down to 0, the test result of control solution is shown on the screen. Compare this result with the range printed on the test strip vial. It should fall within this range.

### Out-of-range results

If test results fall outside the range printed on the test strip vial, check "Operational Problem" section in troubleshooting guide and repeat the test. If you continue to get out-of-range results, it means that the system may not be working properly. Do NOT test your blood. Contact the Customer Care Line at 1-866-469-2632.

## TESTING YOUR BLOOD SUGAR

Be sure to read this section and the test strip package insert found in the test strip box carefully before testing. Make sure you have all items needed to test:


- A. Blood glucose meter      B. Test strip      C. Lancing device  
D. Sterile lancet              E. Clear cap (For AST use)


A.

B.

C.

D.

E.

### Preparing the Lancing Device for Blood Testing

Please refer to the manufacturer's instructions for the lancing device to collect a blood sample.

### WARNING

To reduce the chance of infection:

- The meter and lancing device are for single use. Never share a lancet or the lancing device.
- Always use a new, sterile lancet. Lancets are for single use only.
- Avoid getting hand lotion, oils, dirt, or debris in or on the lancets and the lancing device.
- Wash and dry your hands thoroughly after handling the meter, lancing device, and test strips to prevent infection. For more information, please refer to the "Cleaning and Disinfection" section.
- If the meter is being operated by a second person who is providing testing assistance to the user, the meter and lancing device should be decontaminated prior to use by the second person.


**Sharing the lancing device and lancets may increase the risk of contracting infectious diseases. Lancing device must not be used for more than one person.**

## Test Procedure

### Wash and dry your hands before starting.

#### Step 1. Set the Lancing Device


Screw off the cap of the lancing device. Insert a lancet into the lancet holder and push down firmly until it is fully seated.


Twist the protective disk until it separates from the lancet. Replace the lancet device cap. Turn the cap until it is snug but not too tight.

Twist the adjustable tip in either direction until the number lines up with the arrow: The adjustable tip offers 5 levels of skin penetration:  
1-2 for soft or thin skin, 3 for average skin,  
4-5 for thick or calloused skin.


Slide the ejection/cocking control back until it clicks. If it does not click the device may have been cocked when the lancet was inserted.

#### Blood from sites other than the fingertip

A clear cap, together with the kit, make it easier to get a drop of blood for AST. When you want to obtain blood from sites other than the finger, replace the lancet device cap with the clear cap. Turn the clear cap until it is snug but not too tight, and then slide the ejection/cocking control back until it clicks.


The lancing device is now ready for use. Set aside for later use.


## STEP 2. Insert Test Strip

Fully insert test strip face up with contact bars end first, into the test slot. The Meter will turn on automatically and display the following in sequence:

→ "CH" and "  " → flashing "  " with date and time

## STEP 3. Select the Appropriate Measurement Mode by Pressing M


## STEP 4. Get a Drop of Blood

Select the puncture site either in finger or in other parts (AST). Clean the puncture site with 70% using cotton moistened with 70% alcohol and let it air-dry.

### • Fingertip Testing

Hold the lancing device firmly against the side of your finger. Press the release button. You will hear a click, indicating that the puncture is complete.

After penetration, discard the first drop of blood with a clean tissue paper or cotton. Gently squeeze the punctured area to obtain blood. Be careful NOT to smear the blood sample.

The volume of blood sample must be at least 0.7 microliter in volume (• actual size).

 Correct  Incorrect

- **Sites Other Than Fingertip**


Please refer to the section of “About Alternative Site Testing (AST)” for possible punctured sites.

### PLEASE NOTE

- Choose a different spot each time you test. Repeated punctures at the same spot may cause soreness and calluses.
- Before you decide to begin AST, please consult your health professional.
- It is recommended that you discard the first drop of blood as it might contain tissue fluid and serum, which may affect the test result.


### STEP 5. Apply Blood into the Test Strip

When “” is flashing on the screen, apply your blood to the absorbent hole of the test strip until the confirmation window is **fully covered** with blood. The meter will then begin to count down automatically.

### STEP 6. Obtain an Accurate Result in 7 Seconds

The result of your blood glucose test will appear after the meter counts to 0. This reading will automatically be stored in the meter. Take the strip out from the meter by hand and the meter will show "OFF" then turn off.

### PLEASE NOTE

- Do not push your punctured finger (with blood on it) against the test strip or try to apply a smeared blood sample to the test strip.
- If you do not apply a blood sample to the test strip within 3 minutes, the meter will automatically turn off. You must remove and reinsert the test strip to restart the test procedure.
- The blood should be completely filled the confirmation window before the meter begins to count down. If you find that the confirmation window is not filled with blood when the meter is counting. **NEVER** try to add more blood to the test strip. **Discard the test strip and retest with a new one.**
- If you have trouble filling the confirmation window, contact the Customer Care Line at 1-888-307-8188 for help.


### STEP 7 Remove the used test strip and lancet

Remove the test strip by hand, throw the strip into a disposal container. The meter will turn itself off automatically after the test strip is removed.

**Always use caution when removing the lancet.**

Take the lancet out carefully. Place the disk on a hard surface and push the exposed tip into the protective disk.

**Always follow the instructions in the lancing device insert when removing the lancet.**

### WARNING!

The used lancet and the used test strip may be biohazards. Please discard it carefully according to your healthcare provider's instructions.

## Expected Test Results

---

A Blood glucose meter plays an important role in diabetes control. A long-term study showed that **maintaining normal blood glucose levels** can reduce the risk of diabetes complications by up to 60%.<sup>\*1</sup> The results provided by the FORA V12 No-Coding system can help you and your healthcare professional monitor and adjust your treatment plan to gain better control of your diabetes.

The meter provides you with plasma equivalent results.

Time of day	Normal plasma glucose range for people <b>without</b> diabetes (mg/dL)	Your target range (mg/dL)
Fasting and before meal	Less than 100 mg/dL (5.6 mmol/L)	
2 hours after meals	Less than 140 mg/dL (7.8 mmol/L)	

Reference:

American Diabetes Association (2012). Standards of Medical Care in Diabetes. Diabetes Care, 35 (Supplement 1): S1–S100.

\*1: American Diabetes Association position statement on the Diabetes Control and Complications Trial (1993).

Please work with your doctor to determine a target range that works best for you.

## COMPARING METER AND LABORATORY RESULTS

The meter provides you with plasma equivalent results. The result you obtain from your meter may differ somewhat from your laboratory result due to normal variation. Meter results can be affected by factors and conditions that do not affect laboratory results in the same way (See test strip package insert for typical accuracy and precision data, and for important information on Limitations). To make an accurate comparison between meter and laboratory results, follow the guide-lines below.

### Before you go to the lab:

- Perform a control solution test to make sure that the meter is working properly.
- Fast for at least eight hours before doing comparison tests, if possible.
- Take your meter with you to the lab.

### While at the lab:

Make sure that the samples for both tests (the meter test and the lab test are taken and tested within 15 minutes of each other).

- Wash your hands before obtaining a blood sample.
- Never use your meter with blood that has been collected in a gray-top test tube.
- Use fresh capillary blood only.

You may still have a variation from the result because blood glucose levels can change significantly over short periods of time, especially if you have recently eaten, exercised, taken medication, or experienced stress<sup>\*2</sup>. In addition, if you have eaten recently, the blood glucose level from a finger prick can be up to 70 mg/dL higher than blood drawn from a vein (venous sample) used for a lab test<sup>\*3</sup>. Therefore, it is best to fast for eight hours before doing comparison tests. Factors such as the amount of red blood cells in the blood (a high or low hematocrit) or the loss of body fluid (dehydration) may also cause a meter result to be different from a laboratory result.

### References

\*2: Surwit, R.S., and Feinglos, M.N.: Diabetes Forecast (1988), April, 49-51.

\*3: Sacks, D.B.: "Carbohydrates." Burtis, C.A., and Ashwood, E.R. (ed.), Tietz Textbook of Clinical Chemistry. Philadelphia: W.B. Saunders Company (1994), 959.

## USING THE METER MEMORY


### Viewing Results on the Meter

Your meter stores the 450 most recent blood glucose test results along with their respective dates and times in its memory. It also provides you with 7-, 14-, 21-, 28-, 60- and 90-day averages of your blood glucose test results. You can review these averages test results in memory through these easy steps.

- **Read the Average of Blood Glucose Results:**

When the meter is off, press and release the M button. The screen shows “ ”. Keep pressing the M button for 2-3 seconds, until blinking “ DAY ” appears. Release the M button and then the 7-day average will appear, indicating that you are in the memory mode. If you continue to press the M button, the 14-, 21-, 28-, 60- and 90-day average will appear in order.


The 7-day average is calculated from the blood glucose results obtained during the last 7 days.

It also indicates how many blood glucose tests have been performed within this period, e.g., 14 (14 tests in the last 7 days).


The 14-day average is calculated from the blood glucose results obtained during the last 14 days. It, too, indicates how many blood glucose tests have been performed, e.g., 28 (28 tests in the last 14 days).

The same follows for the 21-, 28-, 60- and 90-day averages.


## • Recall the Stored Test Results

1. When the meter is off, press and release the M button. The screen shows "M". Press the M button again, "01" appears first and then the latest glucose result along with date and time will be shown on the screen.
2. Press the M button once by once to recall the test results stored in the meter consecutively.
3. After the last test result, press the M button again and the meter will be turned off.



## • Exit the Memory Mode

You can also exit memory mode by pressing M button for 5 seconds to turn off the meter.

## PLEASE NOTE

- The control solution results are **NOT** stored in the memory (please refer to **WARNING on the section, CONTROL TEST**, for more information). The list of past results and the result averages are for blood glucose results only.
- When pressing the  button to recall the test results, the meter will first display date and time. If you do not press the  button within 5 seconds, the blinking “” and “” will appear. In the meantime, you can choose to insert a test strip to **start testing your blood** or press the  button again to review **the stored test results**.

- If no button is pressed for 3 minutes, the meter will display “OFF” and turn off automatically.

## Viewing Results on a Personal Computer

---


Results in memory can be transmitted to a personal computer. Health Care System Software and an interface cable are required before installation. The software can be downloaded from the home page of **Fora Care Inc.** To learn more about Health Care System Software or to obtain an Interface Cable separately, please contact the Customer Care Line at 1-888-307-8188 for help.

### Step 1 Install Software

To download the Health Care Software System, please visit the FORA Care Inc. website at <http://www.foracare.com/support.html> and select the "Software & App" to download the "Health Care Software System" software. Follow the online instructions to install the software to your computer.

### Step 2 Connect to Personal Computer

Connect the interface cable to a serial port of your computer. With the meter turned off, connect the Interface Cable to the Data Port of the meter. "PC" will appear on the display, indicating that the meter is ready to transmit data.

### Step 3 Transmit Data

Follow the instructions provided in the software to transmit data. Results transmitted will include date and time. Remove the cable and the meter will automatically turn off.

### PLEASE NOTE


While the meter is connected to the PC, it is unable to perform a blood glucose test.

## MAINTENANCE


### Battery Replacement

---

Your meter comes with two 1.5V AAA alkaline batteries. The meter will alert you when the power is getting low by displaying two different messages:

1. With  symbol appears on the display: the meter is functional and the result remains accurate, but it is time to change the battery.

2. With  symbol, low and E-b symbols on the displays: the battery can not provide enough power for a test. You must change the battery immediately.

#### PLEASE NOTE

- Replacing the batteries within 5 minutes does not affect the settings. Therefore, the settings do not need to be updated when the batteries are replaced.
- Batteries might leak chemicals if unused for a long time. Remove the battery if you are not going to use the device for an extended period (i.e., 3 months or more).
- Do not use the used batteries. Use only the new ones of required size and type.
- Properly dispose the batteries or contact our local customer Care Line for more information.

To replace the battery, make sure that the meter is turned off.

**STEP1**

**STEP1**

Press the buckle on battery cover and lift up to remove cover.

**STEP2**

**STEP2**

Remove the old battery and replace with two 1.5V AAA alkaline batteries.

**STEP3**

**STEP3**

Close the battery cover.

**WARNING**

As with all small batteries, these batteries should be kept away from small children. If they are swallowed, promptly see a doctor for help.

## Caring for Your Meter

---

To avoid the meter and test strips attracting dirt, dust or other contaminants, please wash hands thoroughly with soap and water before and after use.

### **Why the cleaning and disinfection should be performed**

Cleaning and disinfection are different. Cleaning is the process of removing dirt (e.g. food debris, grease, dust), disinfection is the process of killing germs (e.g. bacteria and viruses).

### **When to clean and disinfect the meter**

Clean the meter when you see any dirt on it. You should disinfect the meter at least once a week to prevent infection.

### **NOTE**

If the meter is being operated by a second person who is providing testing assistance to the user, the meter and lancing device should be disinfected prior to use by the second person.

### **How to clean and disinfect the meter**

The meter must be cleaned prior to the disinfection. Use one disinfecting wipe to clean exposed surfaces of the meter thoroughly and remove any visible dirt or blood or any other body fluid with the wipe. Use a second wipe to disinfect the meter. Do NOT use organic solvents to clean the meter.


We recommend using the disinfecting wipes/towelettes below for meter cleaning and disinfection because the active ingredients have been tested to be effective against Hepatitis B Virus (HBV) for **FORA V12 No-Coding blood glucose monitoring system**.

- ▲ Micro-Kill+™ (Micro-Kill Plus™) by Medline (EPA Reg. No. 59894-10-37549)

To obtain disinfecting wipes and other information, please contact Medline at 1-800-MEDLINE (1-888-307-8188) or visit [www.medline.com](http://www.medline.com). You can also purchase at [www.amazon.com](http://www.amazon.com).

## Disinfecting Procedures

1. Take out one disinfecting wipe from the package and squeeze out any excess liquid in order to prevent damage to the meter.

2. Wipe all meter's exterior surface display and buttons. Hold the meter with the test strip slot pointing down and wipe the area around the test slot but be careful not to allow excess liquid to get inside. Keep meter wet with disinfection solution contained in the wipe for a minimum of 2 minutes for Micro-Kill+™ wipes. Follow the instructions on the package label of disinfecting wipe.

3. Remove the wipe. Allow the meter surface to dry completely.

4. Discard the used wipes and never reuse them. Wash your hands thoroughly with soap and water after handling the meter, lancing device and test strips to avoid contamination.

This device has been validated to withstand up to 5000 cleaning and disinfection cycles using the recommended disinfecting wipe/towelette. The tested number of cycles is estimated by 2 cleaning and disinfection cycles per day over 5 years, the expected life of the device. The meter should be replaced after the validated number of cleaning and disinfection cycles or the warranty period, whichever comes first.

Improper system cleaning and disinfection may result in meter malfunction. Stop using the meter if you see any signs of deterioration, for example, LCD display cracks or becomes cloudy, buttons no longer function, or outer casing cracks. If you have any question. Please contact local customer Care Line for a replacement meter if any of the signs of deterioration are noticed.

## NOTE

- Do NOT clean and disinfect the meter while performing tests.
- Please follow the instructions on the package label of Micro-Kill disinfecting wipes for safe use of the wipes
- If the meter is being operated by a second person, the meter and lancing device should be decontaminated prior to use by the second person.
- Do NOT allow cleaning and disinfecting solution to get in the test slot, battery compartment, or strip-ejection button.
- If you do get moisture in the test strip slot, wipe it away with a corner of tissue.
- Always dry the meter thoroughly before using it.
- Do not spray the meter directly with cleaning solutions especially those containing water (i.e. soapy water), as this could cause the solution to enter the case inside and damage the electronic components or circuitry

### **Meter Storage**

- Storage condition: -4°F to 140°F (-20°C to 60°C), below 95% relative humidity.
- Always store or transport the meter in its original storage case.
- Avoid dropping and strong impact.
- Avoid direct sunlight and high humidity.

## Caring for Your Lancing Device

---

The lancing device is reusable. Please disinfect your lancing device regularly.

- To clean the lancing device, use a disinfecting wipe to clean the lancing device.
- Disinfect the cap by placing it in 70% alcohol for 10 minutes and then allow it to air dry.
- Do NOT place the device in a dishwasher or use detergents.

## Caring for Your Test Strips

---

- Storage condition: 39.2°F to 104°F (4°C to 40°C), below 85% relative humidity. Do not freeze.
- Store your test strips in their original vial only. Do not transfer to other container.
- Store test strip packages in a cool and dry place. Keep away from direct sunlight and heat.
- After removing a test strip from the vial, immediately close the vial cap tightly.
- Touch the test strip with clean and dry hands.
- Use each test strip immediately after removing it from the vial.
- Write the opening date on strip vial label when you first opened it. Discard remaining test strips after 3 months.
- Do not use test strips beyond the expiry date. This may cause inaccurate results.
- Do not bend, cut, or alter a test strip in any way.
- Keep the strip vial away from children since the cap and the test strip may be a choking hazard. If swallowed, promptly see a doctor for help.

For further information, please refer to the test strip package insert.

## PROBLEM-SOLVING GUIDE

The following is a summary of various display messages. These messages help to identify certain problems but do not always appear in all cases when a problem has occurred. Improper use may cause an inaccurate result without producing an error message. In the event of a problem, refer to the information under the row "ACTION".



**Never try to disassemble the meter under any circumstances. If you encounter any error messages not listed below or if you have followed the actions recommended below but the problem persists unsolved, please call the Customer Care Line at 1-888-307-8188 for support.**

## Special Message

Special symbols and messages appear together with your test result.

MESSAGE	WHAT IT MEANS	WHAT TO DO
	<p>Your result is below measurement limit, which is less than 20 mg/dL .</p> <p>This indicates hypoglycemia (low blood glucose).</p>	<p>You should seek immediate medical assistance.</p>
	<p>Your result is above measurement limit, which is higher than 600 mg/dL .</p> <p>This indicates severe hyperglycemia (high blood glucose).</p>	<p>You should seek immediate medical assistance.</p>
<b>KETONE?</b>	<p>Your result is equal to or higher than 240 mg/dL .</p> <p>This indicates there is a possibility of ketone accumulation if you are Type 1 diabetes.</p>	<p>Please seek immediate medical assistance.</p>

## Error Message

MESSAGE	WHAT IT MEANS	WHAT TO DO
	Appear when the battery can not provide enough power for a test.	Replace the battery immediately.
	Problem with the meter.	Review the instructions and re-test with a new test strip.
	You need to remove the strip after applying blood to the absorbent hole.	Re-test with a new test strip.
	Appear when environmental temperature is below system operation range.	System operation range is 10°C to 40°C (50°F to 104°F). Repeat the test after the meter and test strip have reached the above temperature.
	Appear when environmental temperature is above system operation range.	
	Appear when inserting a used test strip.	Test with a new test strip.

## Operational Problem

### 1. If the meter does not display a message after inserting a test strip:

PROBABLE CAUSE	WHAT TO DO
Battery exhausted.	Replace the battery.
Battery incorrectly installed or absent.	Check that the battery is correctly installed.
Test strip inserted upside down or incompletely.	Insert the test strip face up with the contact bars end first.
Defective meter.	Please contact the Customer Care Line at 1-888-307-8188.

### 2. If the test does not start after applying the sample:

PROBABLE CAUSE	WHAT TO DO
Insufficient blood sample.	Repeat the test using a new test strip with larger volume of blood.
Defective test strip.	Repeat the test with a new test strip.
Sample applied after automatic shutdown (3 minutes after last user action).	Repeat the test with a new test strip. Apply sample only when  appears on the display.
Defective meter.	Please contact the Customer Care Line at 1-888-307-8188.

### 3. If the control solution test result is out of range:

<b>PROBABLE CAUSE</b>	<b>WHAT TO DO</b>
Error in performing the test.	Read the instructions thoroughly and repeat the test again.
Control solution poorly shaken.	Shake the control solution vigorously and repeat the test again.
Expired or contaminated control solution.	Check the expiration date or the discarded date of the control solution.
Control solution that is too warm or too cold.	Place control solution, meter and test strips in room temperature (68°F to 77°F/20°C to 25°C) before testing.
Test strip deterioration.	Repeat the test with a new test strip.
Meter malfunction.	Please contact the Customer Care Line at 1-888-307-8188.

## SPECIFICATIONS

- Model No.: FORA V12
- Dimension & Weight: 96mm(L) x 45mm(W) x 23mm(H), 76.15g
- Power source: two 1.5V AAA alkaline batteries
- Memory: 450 measurement results with respective dates and times
- Auto electrode inserting detection
- Auto sample loading detection
- Auto reaction time countdown
- Auto shutdown after 3 minutes of idleness
- Operating condition: 50°F to 104°F (10°C to 40°C), below 85% R.H. (noncondensing)
- Meter storage/transportation condition: -4°F to 140°F (-20°C to 60°C), below 95% R.H.
- Strip Storage/Transportation Conditions: 39.2°F to 104°F (4°C to 40°C), below 85% R.H.
- Shelf life: 3 years
- Measurement units: mg/dL
- Measurement range: 20 to 600mg/dL

The device has been tested to meet the electrical and safety requirements of:

IEC/EN 61010-1, IEC/EN 61010-2-101, EN 61326-1, IEC/EN 61326-2-6.

## SUMMARY OF OPERATION

This summary is intended only for quick reference only. It is not a substitute for the owner's manual. Please read the entire manual before you begin testing.

### STEP 1. Insert Test Strip


The meter will automatically turn on and display the followings in sequence:

“CH” and “” → flashing “” with date and time.

### STEP 2. Apply Sample

Touch and hold the drop of blood to the absorbent hole of test strip until confirmation window is completely filled.

### STEP 3. Obtain Accurate Results in 7 Seconds

The meter will begin to count down. The result of your blood glucose test will appear after the meter counts to 0.


### STEP 4. Remove the used test strip and lancet

**Always use caution when removing the lancet.**

Place the disk on a hard surface and push the exposed tip into the protective disk.

Discard the used strip and lancet according to your healthcare provider's instructions.

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