

BINSON'S HOME HEALTH CARE CENTERS

BINSON'S HOME HEALTH CARE CENTERS INCONTINENCE & MEDICAL SUPPLY PROGRAM

Monday through Friday from 8:00 A.M. to 5:30 P.M.

- **Toll Free Number 1-888-217-9610**
 - TTY Relay Service at 711 For Hearing Impaired
- 24-HOUR EMERGENCY SERVICE AVAILABLE

Welcome to Binson's State of Indiana Incontinence and Medical Supply Program. We are committed to providing you with quality incontinence products. All products will need to be ordered by your physician. We also offer an in-house nursing staff for telephone consultation and initial nursing assessments.

Beginning February 1, 2008, Binson's began providing bowel/bladder incontinence supplies by United Parcel System (UPS) or United States Postal System (USPS) through a home delivery contract with OMPP. Your products will be shipped to you monthly unless you report changes.

Our team of customer service representatives is available to answer your questions and assist with your product needs. If you have questions regarding your supplies or did not receive your supplies, please call us at **1-888-217-9610**.

PRODUCT SUPPLIED: UROSTOMY POUCHES & SUPPLIES

DESCRIPTION:

- One or two piece system made of odor-proof, clear or beige plastic with a cloth backing.
- Pouch is attached to a barrier or attaches to a separate barrier that is held to the skin around the stoma with an adhesive wafer.
- Has a spout at the base for emptying.
- Available in a night time form with a tube connection for emptying into a bedside drainage bag.

PRODUCT USE:

Applying/replacing the barrier or barrier/pouch unit

- Gather together soft paper towels or washcloth, warm water, new barrier or barrier/pouch unit, stoma paste (if recommended by you doctor).
- Stand over the toilet or sit far back on the toilet seat. Hold up the end of the pouch and open the spout. Drain contents of the pouch into the toilet bowl.
- Gently remove the stoma/pouch unit or pouch and place in the trash-DO NOT FLUSH.
- To remove the barrier, work slowly from the top of the barrier to the bottom pulling the skin down as you lift off the wafer.

BINSON'S HOME HEALTH CARE CENTERS

Page 2

Patient Education

Urostomy Pouches

- Cleanse the skin around the stoma gently but completely with warm water. Pat to dry or air dry. (Stoma may bleed slightly with friction - this is normal).
- Warm the new wafer between your hand for a minute or two – this will cause better adhesion and longer wearing time.
- Place the wafer around the stoma taking care to avoid wrinkles and securing the wafer with gentle pressure.
- Apply a new pouch (if using a 2 piece unit) after making sure the spout is closed.

EMPTYING THE UROSTOMY POUCH:

- Empty the pouch into the toilet when approximately one thirds to one half full.
- Stand over the toilet or sit as far back as possible on the toilet seat.
- Hold up the end of the pouch and open the spout.
- Drain into the toilet.
- Wipe off the spout thoroughly before closing.

SKIN CARE:

- Inspect the stoma and the surrounding skin with every wafer wafer/pouch change.
- Wash the area with warm water. Cleansers without oils, glycerin or deodorants may be used.
- Stoma should be a healthy, rosy red color.
- Report any rashes or skin problems to your doctor or ostomy nurse immediately.