

BINSON'S HOME HEALTH CARE CENTERS

BINSON'S HOME HEALTH CARE CENTERS INCONTINENCE MEDICAL SUPPLY PROGRAM

Monday through Friday from 8:00 A.M. to 5:30 P.M. EST

- **Toll Free Number 1-888-217-9610**
 - TTY Relay Service at 711 For Hearing Impaired
- 24-HOUR EMERGENCY SERVICE AVAILABLE

Welcome to Binson's State of Indiana Incontinence and Medical Supply Program. We are committed to providing you with quality incontinence products. All products will need to be ordered by your physician. We also offer an in-house nursing staff for telephone consultation and initial nursing assessments.

Beginning February 1, 2008, Binson's began providing bowel/bladder incontinence supplies by FedEx, United Parcel Service (UPS) or United States Postal Service (USPS) through a home delivery contract with Indiana Medicaid. Monthly, your products will be shipped to you monthly unless you report changes.

Our helpful team of customer service representatives is available to answer your questions and assist with your product needs. If you have a question regarding your supplies or you did not receive your supplies, please call us at **1-888-217-9610**.

PRODUCT SUPPLIED: **OSTOMY ACCESSORIES**

DESCRIPTION:

- Filters, deodorants, clamps, belts, absorbent thickener, tape, paste, irrigation supplies or cleaners intended for cleaning, deodorizing, anchoring/adhering or managing an ostomy.

PRODUCT USE:

- **Filters** are used to absorb odor and are often made of charcoal. Has adhesive strip to adhere to pouches. Follow manufacturer's package instructions.
- **Deodorants** can be drops or tablets that are put in the pouch to eliminate odor. Follow manufacturer's package directions.
- **Clamps** are used to secure the open end of a drainable pouch. Most commonly made of plastic, are curved to conform to body with a simple click-type closure.
- **Belts** are used to secure the position of the barrier/flange to the abdomen wall and strengthen the seal between the barrier and skin. Made of cotton, rayon and a small amount of latex for flexibility and are adjustable. Hand wash when soiled and line dry. Follow manufacturer's package instructions.
- **Absorbent thickener** is usually a cornstarch-based powder sprinkled into the pouch to increase the thickness of drainage. Follow manufacturer's package instructions.
- **Tape** can be non-waterproof or waterproof and is used to secure a barrier.
- **Paste** can be in a tube or moldable strips. It is used between the body and the barrier to fill in gaps and create a flat surface to accept the barrier. Should be applied to clean dry skin. May contain alcohol - Do not apply if the skin around the stoma is reddened or broken. Follow manufacturer's package directions.
- **Irrigation supplies** are used to irrigate the colon and loosen stool or increase its fluid content. Collection bags, tubing, stoma cone adapter and faceplate are all disposable. Follow manufacturer's package directions.
- **Cleaners** are liquids used to clean and freshen a reusable pouch or appliance of urine/fecal matter and deodorize at the same time. Follow manufacturer's package directions.

SKIN CARE:

- View the skin daily. Report any redness, rash or breaks in the skin to your physician.

BINSON'S HOME HEALTH CARE CENTERS

- Clean the skin daily and after soiling. Use warm water with mild cleanser.
- Pat the skin dry or air dry before applying a new barrier.