



# Medline Industries, LP

## **IMMEDIATE ACTION REQUIRED**

### Medline Red Rubber Intermittent Catheters(8Fr-22 Fr) and Medline Red Rubber Intermittent Catheter Trays

05/04/2026

Dear Valued Medline Customer:

Medline Industries, LP is issuing a recall for specific item(s) and lot(s) of Medline Red Rubber Intermittent Catheters (8 Fr-22 Fr) and Medline Red Rubber Intermittent Catheter Trays. During an internal review, Medline identified representative samples did not meet applicable performance specifications. Identified issues include irregularities in catheter surface finish, instances where the catheter connector detached prematurely from the catheter funnel, and limited cases where the catheter outer diameter did not meet specification. Irregular or inconsistent surfaces of the catheter may result in increased friction, which may lead to added discomfort and minor disruptions to the intermittent catheter insertion. Separately, catheter outer diameters outside of specification may also contribute to added discomfort and minor disruption during insertion. With respect to the catheter connector, premature catheter disconnection and leakage of urine may disrupt a medical procedure and potentially prompt discontinuation and re-catheterization. All identified hazards are not expected to have more severe or long-term health consequences with use. Inherent risks associated with catheter use, such as urinary tract infection and urinary retention, are not expected to be impacted or elevated by the catheter issues identified. There have been no reports of illness or injury associated with these issues. All affected Medline Sterile Red Rubber Intermittent Catheters and Red Rubber Intermittent Catheter Trays must be destroyed. Be advised that there may be item number variations based on the unit of measure ordered by your facility, indicated by an H appended to the item number(s) included with this communication. Note, there is no product difference between variations in units of measure and physical labeling will reflect the base item number (i.e. DYND13508). Please refer to the enclosed list for all affected item number(s) and affected lot number(s).

Adverse reactions or quality problems experienced with the use of this product may be reported to the FDA's MedWatch Adverse Event Reporting program either online, by regular mail or by fax and to Medline Industries, LP.

#### **REQUIRED ACTION:**

1. Immediately check your stock for the affected item number(s) and the affected lot number(s) listed on the enclosed response form. Quarantine all affected product. Upon completion of the below actions, please destroy affected product.
2. Please complete and return the enclosed response form listing the quantity of affected product you have in inventory. Even if you do not have any affected product in inventory, please complete and return the response form.
3. Your account will receive credit once the response form is submitted, if applicable.
4. If you are a distributor, or have resold or transferred this product to another company or individual, you are required by FDA regulations to notify them of this recall communication. Please request that customers destroy any affected product. You should include your customers' quantities on your response form.

If you have any questions, contact the Recall Department at 866-359-1704 or [recalls@medline.com](mailto:recalls@medline.com).

Please accept our sincere apologies for any inconvenience this may have caused. We, like you, place the health and safety of your patients first and foremost.

Sincerely,

**Karin Johnson**

Karin Johnson  
Product Recall Specialist Sr.

Recall Ref#: R-26-080

Confidential - For intended recipients only. Do not Distribute.

Medline Industries, LP. Three Lakes Drive Toll Free: 866.359.1704 Email:  
Northfield, IL 60093-2753 Website: medline.com [Recalls@Medline.com](mailto:Recalls@Medline.com)  
Fax: 866.767.1290



### Item and Lot Numbers

<b>Item Number</b>	<b>Lot Numbers</b>
DYNC1816	All lot Numbers
DYND10350	All lot Numbers
DYND10407	All lot Numbers
DYND13508	All lot Numbers
DYND13510	All lot Numbers
DYND13512	All lot Numbers
DYND13514	All lot Numbers
DYND13516	All lot Numbers
DYND13518	All lot Numbers
DYND13520	All lot Numbers
DYND13522	All lot Numbers
DYND13614	All lot Numbers
DYND13616	All lot Numbers
DYND18350	All lot Numbers
URO4RR14C	All lot Numbers
URO4RR16C	All lot Numbers
DYND13515	All lot Numbers

Confidential - For intended recipients only. Do not Distribute.

**Medline Industries, LP.** Three Lakes Drive      Toll Free: 866.359.1704      Email: [Recalls@Medline.com](mailto:Recalls@Medline.com)  
Northfield, IL 60093-2753      Website: [medline.com](http://medline.com)      Fax: 866.767.1290